# BIOQUANT

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
## Volume Measurements

|            |                                       |
|------------|---------------------------------------|
| TV         | Tissue Volume (mm <sup>2</sup> )      |
| BV         | Bone Volume (mm <sup>2</sup> )        |
| Md.V       | Mineralized Volume (mm <sup>2</sup> ) |
| OV         | Osteoid Volume (mm <sup>2</sup> )     |
| BV/TV      | Bone/Tissue Volume (%)                |
| Md.V/BV    | Mineralized/Bone Volume (%)           |
| Md.V/TV    | Mineralized/Tissue Volume (%)         |
| OV/BV      | Osteoid/Bone Volume (%)               |
| OV/TV      | Osteoid/Tissue Volume (%)             |
| Additional | Unlimited Custom Indices              |

## Surface Measurements

|            |                                       |
|------------|---------------------------------------|
| BS         | Bone Surface (mm)                     |
| OS         | Osteoid Surface (mm)                  |
| Ob.S       | Osteoblast Surface (mm)               |
| ES         | Erosion Surface (mm)                  |
| Oc.S       | Osteoclast Surface (mm)               |
| QS         | Quiescent Surface (mm)                |
| sLS        | Single Labeled Surface (mm)           |
| dLS        | Double Labeled Surface (mm)           |
| BS/BV      | Bone Surface / BV (mm <sup>-1</sup> ) |
| BS/TV      | Bone Surface / TV (mm <sup>-1</sup> ) |
| OS/BS      | Osteoid/Bone Surface (%)              |
| ES/BS      | Erosion/Bone Surface (%)              |
| QS/BS      | Quiescent/Bone Surface (%)            |
| Ob.S/BS    | Osteoblast/Bone Surface (%)           |
| Oc.S/BS    | Osteoclast/Bone Surface (%)           |
| Additional | Unlimited Custom Indices              |

## Thickness Measurements

|            |                            |
|------------|----------------------------|
| O.Th       | Osteoid Thickness (mcm)    |
| W.Th       | Wall Thickness (mcm)       |
| Ir.L.Th    | Interlabel Thickness (mcm) |
| Additional | Unlimited Custom Indices   |

## Number Measurements

|            |                                      |
|------------|--------------------------------------|
| N.Ob       | Osteoblast Number                    |
| N.Oc       | Osteoclast Number                    |
| N.Ob/OS    | Osteoclast Numb/Osteoid Surf (Ob/mm) |
| N.Oc/ES    | Osteoclast Numb/Erosion Surf (Oc/mm) |
| Additional | Unlimited Custom Indices             |

## Structural Indices

|            |   |
|------------|---|
| Tb.Th      | Trabecular Thickness (Plate Model) (mcm)            |
| Tb.N       | Trabecular Number (Plate Model)(mcm <sup>-1</sup> ) |
| Tb.Sp      | Trabecular Separation (Plate Model) (mcm)           |
| Additional | Unlimited Custom Indices                            |

## Kinetic Indices

|            |                                    |
|------------|------------------------------------|
| MS         | Mineralizing Surface (mm)          |
| MS/BS      | Mineralizing/Bone Surface (%)      |
| MAR        | Mineral apposition rate (mcm/day)  |
| Aj.AR      | Adjusted apposition rate (mcm/day) |
| BFR/BS     | Bone formation rate (mcm/day)      |
| Mlt        | Mineralization lag time (days)     |
| Omt        | Osteoid maturation time (days)     |
| FP         | Formation period (days)            |
| Rs.P       | Resorption period (days)           |
| Rv.P       | Reversal period (days)             |
| Rm.P       | Remodeling period (days)           |
| QP         | Quiescent period (days)            |
| Tt.P       | Total period (days)                |
| Ac.f       | Activation frequency (/year)       |
| Additional | Unlimited Custom Indices           |

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

D3 Tissue Volume, microns squared  
D4 Bone Volume, microns squared  
D5 Tumor Volume, microns squared  
D6 Osteoclast Number, number (count)  
L1 Bone Surface, microns  
L2 Osteoclast Surface, microns  
L3 Tumor Surface, microns

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BV/TV: Bone to Tissue Volume Ratio, percent  
BV/Tu.V: Bone to Tumor Volume Ratio, percent  
BS/BV: Bone Surface to Bone Volume, 1/mm  
Oc.S/BS: Osteoclast Surface to Bone Surface, 1/mm  
Oc.N/Oc.S: Osteoclast Number to Osteoclast Surface, #/mm  
Oc.N/BS: Osteoclast Number to Bone Surface, #/mm  
Tu.S/BS: Tumor Surface to Bone Surface, 1/mm

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PCNA labeled Section.  
D1 Topo  
D4 Proliferation Zone  
D5 PCNA+ Chondrocytes, microns squared  
D6 PCNA+ Chondrocytes, microns squared

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
Pr.Z.Ar: Total Proliferation Zone Area,  $\text{mm}^2$   
PCNA+ Ch.N: PCNA+ Chondrocyte Number, number  
PCNA+ Ch.N: PCNA+ Chondrocyte Number, number  
T.Ch.N: Total Chondrocyte Proliferation, number  
PI: Proliferation Index  
+Ch.N/Pr.Z.AR: PNCNA+Chondrocyte Number/  
Proliferation Zone, number/ $\text{mm}^2$

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



## Fluorescent Indices

|                              |                            |
|------------------------------|----------------------------|
| (Ec.M.Pm)                    | Ec. Mineralizing Pm        |
| (Ec.Ir.L.Wi)                 | Avg Interlabel Wi          |
| (Ec.MAR)                     | Mineral Apposition Rate    |
| (Ps.Ir.L.Wi)                 | Avg Interlabel Wi          |
| (Ps.MAR)                     | Mineral Apposition Ra      |
| (Ps.sL.Pm)                   | Periosteal Single Label Pm |
| (Ps.dL.Pm)                   | Periosteal Double Label Pm |
| (Ps.M.Pm)                    | Periosteal Mineralizing Pm |
| Unlimited Additional Indices |                            |

## Formation Indices

|                              |                      |
|------------------------------|----------------------|
| (Ac.F)                       | Activation Frequency |
| (BFR)                        | Bone Formation Rate  |
| Unlimited Additional Indices |                      |

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Edit Arrays Image Measure Topo Options BQ Scan Stereo Topographer Help  
 Welcome  
 come to BIOQUANT OSTEO! This is the Guide  
 ei. It will help introduce you to BIOQUANT. To  
 tinue, click the Next button.

Type: Image File  
 Open Single  
 Open Sequence  
 Save Image  
 Calibration  
 Clean  
 Redraw

New Set  
 Quick Set  
 Modify  
 Calculate  
 Open List  
 Save List

Available Arrays  
 D13 Osteoid Volume (0)  
 D14 Adipocyte Volume (0)  
 D15 Osteoblast Number (0)  
 D16 Osteoclast Number (0)  
 L1 Osteoid Surface (0)  
 L2 Osteoclast Surface (0)  
 L4 Osteoblast Surface (0)

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Save to Topo Array:  
 Mag: X20  
 Topo Array: D4  
 Z Offset: 0  
 Color (C): [Color Picker]


Type: Topo ROI  
 Define Assign  
 Landmark Move

Refresh Show Transparent  
 Assign  
 Add Undo

Threshold  
 R G B  
 [Color Bars]  
 R[136, 43], G[193, 90], B[204, 119]

Type: General  
 Select  
 All Objects  
 Assign Measure Preview Seg Assign


Available Filters  
 Void Filter  
 Outline Filter  
 Measurement

Calculations  
 File Data Update  
 A B C D E

|                |       |
|----------------|-------|
| P1 TV          | 5.33  |
| P2 BV          | 0.84  |
| P3 BV/TV       | 0.16  |
| P4 BS          | 23.54 |
| P5 BS/BV       | 28.12 |
| P6 Tb.Dm (Red) | 0.14  |
| P7 Tb.N (Red)  | 3.14  |
| P8 Tb.Sp (Red) | 0.18  |
| P9 Adv         | 0.00  |
| P10 Adv/TV     | 0.00  |
| P11 OV         | 0.00  |
| P12 OV/BV      | 0.00  |
| P13 OS         | 0.00  |
| P14 OS/BS      | 0.00  |
| P15 O.VI       | 0.00  |
| P16 Oc.S       | 0.00  |
| P17 Oc.S/BS    | 0.00  |
| P18 N.Oc       | 0.00  |
| P19 N.Oc/BS    | 0.00  |
| P20 N.Oc       | 0.00  |
| P21 N.Oc/BS    | 0.00  |

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邮箱：carlyang@tr-baast.com 网站：www.tr-baast.com

---

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TONGREN IMPORT AND EXPORT TRADING ( HONG KONG ) LIMITED

Add:Room A, 7/F, China Overseas Building 139 Hennessy Road,

Wanchai,Hong Kong

Tel: ( 852 ) 28511477 Fax: ( 852 ) 28511480

邮箱：yanglq811@163.com

---

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BIOQUANT Image Analysis Corporation

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